

# Improved Productivity Through Continuous Workflow

Automatic Load and Unload System for  
Amsco® Evolution™ Medium Steam Sterilizers



The Automatic Load and Unload System (ALUS) from STERIS can help save time and improve efficiencies so staff can focus on other tasks.

- Automatically loads and unloads carts when cycle is complete
- Starts cycle automatically with optional bar code reader

## Continuous Workflow

The next load can be set up while another load is being sterilized. After the cycle has finished, the system automatically moves the sterile load onto the unload table, where it can stay until staff is available to transport it for storage or distribution.

## High-Volume Capacity

The system's heavy-duty construction can handle loads of up to 500 pounds per load (in the largest chamber models).

## Ergonomic Design

Use the transfer carriage to safely move the loading cart onto the load table, and remove it from the unload table, using minimal effort – which helps to reduce the risk of injuries.

## Intuitive Controls

Push-button controls are located on the side of the each table and the loading operation is fully integrated with the sterilizer controls.

Use the sterilizer display screen to check the load status.



## Choose from Three Options

- ALUS without a bar code reader
- ALUS with a bar code reader
- Manually load and automatically unload without a bar code reader

The optional bar code reader determines the appropriate preprogrammed cycle, which starts automatically once the load is inside the Evolution Medium Steam Sterilizer.

## Specifications

- 26" x 26" or 26" x 37-1/2" chamber
- Automated chamber doors



Specially designed cart helps facilitate the transfer process.

### Mixed Sources

Product group from well-managed forests, controlled sources and recycled wood or fiber



Printed with Soy based inks and water based coatings on FSC/Green Seal-certified paper.

Document # M3444EN.2010-06, Rev. A  
Printed 06/2010, 2000

©2010 STERIS Corporation.  
All rights reserved. Printed in USA.

 **STERIS**

STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 · USA  
440-354-2600 · 800-548-4873  
[www.steris.com](http://www.steris.com)