

## AMSCO® C Series

Certified Pre-Owned  
Small Steam Sterilizers

**CERTIFIED**  
STERIS multi-point checklist

**COMPLIANT**  
AAMI-ST8:2008: ST79: 2010

**CAPACITY**  
Processes two 25 lb trays

AMSCO® C Series  
STERIS



# AMSCO C is the only refurbished steam sterilizer upgraded to meet current AAMI standards.



## Certified pre-owned

All STERIS Certified pre-owned equipment:

- Is reviewed and refurbished using STERIS original quality assurance guidelines, multi-point checklist and original STERIS parts, as required
- Passes the same factory test procedures as new STERIS equipment
- Comes with a warranty

AMSCO C Steam Sterilizers are upgraded to current compliance, safety and technology standards.

## Compliant

The AMSCO C Steam Sterilizer includes a new control system that meets the most current AAMI standards for immediate-use sterilization and has the capacity to process today's heavier instrument sets.

- AAMI-ST8: 2008
- AAMI-ST79:2010
- Meets UL/CSA/IEC safety standards

## Multiple configurations

- Single or double door
- Cabinet enclosed/freestanding or recessed
- House steam
- Integral generator with flush and drain

## Accessories

- Rack and two shelves for 20" units (included with 16" configuration)
- Loading car/ transfer carriage (20" only)

**Capacity** The AMSCO C Steam Sterilizer processes two 25 lb trays, which means it's more productive.\*

Model	Size	Type	Load	Max Weight
AMSCO C Series 16C	16x16x26"	Prevac	2 trays/4 fabric	50 lbs
AMSCO C Series 20C	20x20x38"	Prevac	2 trays/6 fabric	50 lbs
AMSCO C Series 16CS	16x16x26"	SFPP	2 trays/4 fabric	50 lbs
AMSCO C Series 20CS	20x20x38"	SFPP	2 trays/6 fabric	50 lbs



Visit <http://www.steris.com/products/pre-owned/steam-sterilizers> or contact your STERIS representative to learn more.



STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 · USA  
440-354-2600 · 800-548-4873  
[www.steris.com](http://www.steris.com)

\*Compared to previous generation (AMSCO® Century® Steam Sterilizer)

Document # M4370EN.2016.04, Rev. C  
©2016 STERIS Corporation.  
All rights reserved. Printed in USA.