

More Productive. More Intuitive. More Green.

Vision® 1300 Series Cart and Utensil Washer/Disinfector



The Vision 1321 and 1327 cart and utensil washer/disinfectors deliver maximum throughput and unparalleled performance because they are:

- More Productive – cycle times as fast as 10 minutes
- More Intuitive – effortless operation and reduced downtime
- More Green – utility-saving technology and low cost of ownership



More Productive


Revolutionary Performance

The unique HexaClean spray pattern technology and the enhanced drying system provide more effective cleaning from start to finish:

- high-impingement HexaClean spray pattern provides continuous full-contact washing throughout the entire cycle
- 10-minute cycle times process more carts in less time (50% more productive than competition)*
- provides thermal disinfection
- enhanced drying system allows faster drying time

The Vision 1321 and 1327 cart and utensil washers are manufactured in compliance with the latest North American and global standards.

Case carts processed per 8 hour day*



HexaClean spray in action
Powerful cleaning, proven results

Instrument Processing Option

The optional instrument cycle package provides an alternative for reprocessing instrument trays and loaner sets. This package includes:

- instrument rack with integrated spray arms (validated and compliant specifically with Vision 1300 series washers)
- enzymatic solution pump system
- integrated flow meters to ensure dosing accuracy
- instrument cycle pre-programmed into control
- intermediate level disinfection



Integrated manifold recognizes that instrument rack has been loaded.

*Based on published data of 15 minute cycle times.




More Intuitive

Easy-to-use Control Panel

- single-touch icon buttons for easy cycle selection
- cycle phase and time is visible across the room
- provides cycle parameter verification
- convenient record keeping capabilities
 - USB port – download cycle records at the touch of a button
 - integrated cycle printer on clean side

Smart Technology

- automatic sliding glass doors are easy to operate
- hydraulic seal maintains heat inside the chamber
- automatic floor tilt helps provide more efficient drying
- safety cable on the interior of the chamber
- streamlined, easy-access mechanical core allows for quick maintenance to help maximize uptime
- ProConnect Remote Monitoring technology watches every cycle to make sure your equipment operates consistently.



USB port for downloading records



Automatic single-touch sliding glass door




More Green

The Vision cart washer has been designed with water-and energy-efficient features to provide the lowest cost of ownership (up to 68% savings** compared to the competition).

The Vision 1300 Series:

- provides the lowest water consumption per fill of any washer on the market
- recycles 70% of water each cycle (excludes instrument cycle)
- Prolystica Cleaning Chemistries are formulated using biodegradable ingredients and use recyclable packaging
- 10 minute cycles allow for low utility consumption
- full glass door is made with 20% recycled glass

Gallons of water used per fill



**Calculations based on average national utility costs.

Racks for Every Need

Newly designed racks for the Vision 1300 series washers offer:

- maximum capacity to meet your workflow needs
- easy to load racks that optimize cleaning efficacy
- ergonomic, cool touch handles for safety and efficiency
- convenient removable basket (on container rack) makes loading retention filters easy



Container Rack



Utensil Rack



Cool touch handle

Maximize Cleaning Performance

Efficient cleaning starts with the right chemistries. The lightweight, ergonomic Prolystica Cleaning Chemistries:

- are specifically formulated to maximize the performance of Vision washers
- use less chemistries than non-concentrated formulations
- are available in 2x or 10x concentrations
- have biodegradable ingredients



Value Beyond the Product

ProConnect™ Remote Monitoring System

- 24/7 equipment monitoring
- proactive equipment alerts and issue resolution
- priority technical phone support



To see a video of the cart washer in action, scan this QR code or visit steris.com/visioncartwasher

We're doing our part by printing on recycled paper.

Document #M3533EN.2011-08, Rev. A
Printed 10/2012, 1000



©2011-2012 STERIS Corporation.
All rights reserved. Printed in USA.

 **STERIS**

STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
440-354-2600 · 800-548-4873
www.steris.com